



# Radiation Safety for Staff in Fluoroscopy Suites


Beth Schueler  
Mayo Clinic  
Rochester, Minnesota

## Learning Objectives

- For staff performing fluoroscopically-guided interventional procedures:
  - What are typical radiation exposure levels?
  - How should the radiation exposure to staff be monitored?
  - What type of radiation safety education is needed?
  - Is there anything new and novel available that can help reduce staff exposure levels?

## Resources

- Joint SIR / CIRSE Guideline for Occupational Radiation Protection in IR (Miller et al, 2010)
- NCRP Report No. 168 Radiation Dose Management for Fluoroscopically-guided Interventional Medical Procedures, 2011



## Operator Exposure During Fluoroscopy Procedures


Procedure Type	Mean Dose per Procedure ( $\mu\text{Sv}$ )			Reference
	Neck	Lens	Hand	
Mixed general IR	30-325	300	270-400	NCRP 168
ERCP	450	550	640	NCRP 168
Endovascular surgery	300		400	NCRP 168
Percutaneous coronary intervention	10-130	10-170	30-350	Kim et al, 2008
Cardiac ablation	8-200	50-320	40-230	Kim et al, 2008

## Typical Operator Exposure Levels

- Annual doses for a workload of 1000 procedures
  - Neck: 10-450 mSv
  - Lens of the eye: 10-550 mSv
  - Hand: 30-640 mSv
- Survey of interventional radiologists with a mixed workload (Marx et al, 1992)
  - Mean annual dose (dosimeter on chest over protective apron): 49 mSv (range: 3-115 mSv)

5

## Mayo Clinic Rochester



6

## Operator Exposure Levels

- Typical doses for operators performing fluoroscopically-guided interventional procedures are high
  - May exceed annual dose limits
    - Lens of the eye: 150 mSv
    - Hands: 500 mSv
  - Values generally well exceed those for other healthcare workers
    - Annual US healthcare workers (NCRP Report No. 160): 81% with recordable dose are  $< 1$  mSv

7

## Personnel Dose Monitoring

- Due to the potential for high occupational doses, appropriate monitoring is critical
- Monitoring considerations:
  - 1-dosimeter or 2-dosimeter monitoring?
  - Dose calculation method?
  - When should dose readings be investigated?

8

## Personnel Dose Monitoring



1-Dosimeter Method:  
At neck, outside apron



2-Dosimeter Method:  
At neck, outside apron and  
at waist or chest, under apron

9

## Personnel Dose Monitoring

- 2-Dosimeter method recommended (NCRP Report No. 168)
  - Provides an indication of apron attenuation
  - Allows for better estimate of operator effective dose
  - Dosimeters should be clearly labeled to avoid mix-ups

10

## Dose Assignment with Protective Aprons

- Various calculation models exist and state regulations vary
- CRCPD SSRs (Webster's method)
  - $H_E = 1.5 \times \text{under-lead reading} + 0.04 \times \text{neck over-lead reading}$ 
    - $H_E$ , effective dose equivalent
- NCRP Report No. 122 recommendation
  - $E = 0.5 \times \text{under-lead reading} + 0.025 \times \text{neck over-lead reading}$ 
    - $E$ , effective dose

11

## Operator Exposure During Fluoroscopy Procedures

Procedure Type	Mean Annual Effective Dose* (mSv)	Reference
Mixed general IR	2 - 15	NCRP 168
ERCP	21	NCRP 168
Endovascular surgery	23	NCRP 168
Percutaneous coronary intervention	0.2 - 9	Kim et al, 2008
Cardiac ablation	0.2 - 10	Kim et al, 2008

\* assuming workload of 1000 procedures per year

12

## Annual Dose Limits

- NCRP Report No. 116 recommendation:
  - Effective dose: 50 mSv
- ICRP Publication 60 recommendation:
  - Effective dose: 20 mSv, averaged over 5 years, not to exceed 50 mSv in a single year

13


## Lens Exposure

- Annual dose limit recommendations (NCRP Report No. 116):
  - Lens of the eye: 150 mSv
- Measurements for an annual workload of 1000 procedures,
  - Lens of the eye (unprotected): 10-550 mSv

14

## Lens Dose Estimation

- Exposure level at the eye is typically somewhat lower than at the neck
  - Lens:Neck dose ratio varies with C-arm angulation
- Neck dosimeter reading provides conservative estimate (NCRP No. 168)



15

## Leaded Eyewear

- Typical lead equivalent thickness of radiation protective eyewear is 0.75 mm
  - 98% attenuation
- Actual lens dose is higher due to
  - Exposure from the side and from below
  - Backscatter from head



16



## Leaded Eyewear

- Traditional style

- 0.75 mm lead equivalent lenses
- 120 g
- 28 cm<sup>2</sup> surface area



## Leaded Eyewear

- Sport-wrap style

- 0.75 mm lead equivalent lenses
- 59 g
- 16 cm<sup>2</sup> surface area




## Leaded Eyewear Attenuation

Leaded Eyewear Style	Attenuation Factor		
	0° Angle	45° Angle	90° Angle
Traditional	10	5.2	4.0
Sport-wrap	8.3	4.5	1.4

## New Guideline on Lens Exposure

- ICRP issued a new recommendation (ICRP, 2011)
  - Lower threshold for cataract formation: 0.5 Gy (previous threshold 2-5 Gy)
  - Lower occupational eye dose limit: 20 mSv/yr averaged over 5 years with no year > 50 mSv



## Radiation-Induced Cataract

- Problems with earlier studies:
  - Short follow-up period – latency period is longer for low doses
  - Insufficient sensitivity to detect early lens changes
  - Few subjects with doses below a few gray
- Significant studies:
  - Chernobyl nuclear reactor accident cleanup workers (Worgul et al, 2007)
  - US radiologic technologists (Chodick et al, 2008)

21

## Hand Dose

- Ring dosimeters recommended if hand dose  $> 50$  mSv in a year (NCRP Report No. 168)
- Monitor for a trial period of several months for new staff and new procedure types
- Wear with sensitive area toward exposure source
  - Inward for under-table x-ray tube configurations



22

## Investigation of Dose Readings

- Recommended investigation trigger level (WHO, 2000):
  - E, effective dose > 0.5 mSv/month
  - Lens dose > 5 mSv/month
  - Hand dose > 15 mSv/month
- Verify validity of measurement
- Look for changes in procedure volume, procedure type, equipment, ...

23

## Investigation of Dose Readings

- It is common for personnel who may receive a high occupational dose to not wear their dosimeters to avoid investigations
  - 43% of surveyed interventional radiologists indicate they rarely or never wear monitoring dosimeters (Marx et al, 1992)
- Dosimeter readings that are lower than expected for a specific work assignment should also be investigated (NCRP Report No. 168)

24

## Learning Objectives

- For staff performing fluoroscopically-guided interventional procedures:
  - What are typical radiation exposure levels?
  - How should the radiation exposure to staff be monitored?
  - What type of radiation safety education is needed?
  - Is there anything new and novel available that can help reduce staff exposure levels?

25

## Radiation Safety Education

- Credentials and privileges for fluoroscopy operators are needed
  - Help for developing a program is coming
  - AAPM TG 124 (Chair – Mary Moore)
    - “A Guide for Establishing a Credentialing and Privileging Program for Users of Fluoroscopic Equipment in Healthcare Organizations” is in the works
    - Will include suggestions to encourage your facility to approve a program, didactic content, evaluating competency and information resources for teaching

26


## Occupational Radiation Safety Resources

### ■ IAEA RPOP web pages

- <https://rpop.iaea.org/RPOP/RPoP/Content/InformationFor/HealthProfessionals/index.htm>

### ■ IAEA slide presentation series

- [https://rpop.iaea.org/RPOP/RPoP/Content/AdditionalResources/Training/1\\_TrainingMaterial/Radiology.htm](https://rpop.iaea.org/RPOP/RPoP/Content/AdditionalResources/Training/1_TrainingMaterial/Radiology.htm)



27


## Occupational Radiation Safety Resources

### ■ Image Gently pediatric IR presentation

- <http://www.pedrad.org/associations/5364/ig/>

### ■ RSNA/AAPM Online physics modules

- <http://physics.rsna.org/default.asp>
- “Radiation Safety and Dose in Interventional Radiology”



28

## Learning Objectives

- For staff performing fluoroscopically-guided interventional procedures:
  - What are typical radiation exposure levels?
  - How should the radiation exposure to staff be monitored?
  - What type of radiation safety education is needed?
  - Is there anything new and novel available that can help reduce staff exposure levels?

29

## Operator Shielding

- Overhead and table shields can be very effective
- But may be cumbersome for certain procedures:
  - C-arm angulation
  - Biliary or transjugular access




30



## Orthopedic Complications from Lead Apron Use

- Back pain was reported by 50-75% of interventional physicians surveyed (Klein et al, 2009)
  - Compare to typical incidence of 27% in US adults
  - 25-30% reported that back problems had limited their work
- Options for relief
  - Lightweight aprons
  - Vest/kilt design



31

## Radiation Protective Cabins

- ZeroGravity
  - 1.25 mm lead apron and 0.5 mm lead-equivalent face shield



Marichal et al, 2011

32



## Radiation Protective Cabins

### ■ CATHPAX

- 2 mm lead walls and lead-equivalent windows



Dragusin et al, 2007

33

## Real-time Personnel Dose Monitoring

### ■ DoseAware

- Displays cumulative dose and dose rate on a monitor
- Can be networked between multiple procedure rooms
- Allows for real-time feedback to avoid high scatter conditions and implement radiation reduction techniques



Sanchez et al, 2010

34

## Heavy Metal Protective Patient Drapes

### ■ RADPAD

- Gel pad with tungsten-antimony
- Sterile, dispose after procedure
- 12× eye dose reduction
- 29× hand dose reduction



Dromi et al, 2006

35

## References

- Chodick G, Bekiroglu N, Hauptmann M, et al. Risk of cataract after exposure to low doses of ionizing radiation: A 20-year prospective cohort study among US radiologic technologists, *Am. J. Epidemiol.* 2008; 168(6), 620–631.
- Conference of Radiation Control Program Directors (CRCPD) Suggested State Regulations, <http://www.crcpd.org/ssrcr.aspx>
- Dragusin O, Weerasooriya R, Jais P, et al. Evaluation of a radiation protection cabin for invasive electrophysiological procedures, *Eur. Heart J* 2007; 28(2), 183–189.
- Dromi S, Wood BJ, Oberoi J, Neeman Z. Heavy metal pad shielding during fluoroscopic interventions. *J Vasc Interv Radiol* 2006; 17:1201–1206.
- International Commission on Radiological Protection (ICRP), ICRP Publication 60, Ann. ICRP 1990, 21(1-3) (Elsevier, New York).
- International Commission on Radiological Protection (ICRP), Statement on tissue reactions. 2011. <http://www.icrp.org/page.asp?id=123>
- Kim KP, Miller DL, Balter S, et al. Occupational radiation doses to operators performing cardiac catheterization procedures. *Health Phys* 2008; 94:211–227.
- Klein LW, Miller DL, Balter S, et al. Occupational health hazards in the interventional laboratory: time for a safer environment. *J Vasc Interv Radiol* 2009; 20:147–152.
- Marichal DA, Anwar T, Kirsch D, et al. Comparison of a suspended radiation protection system versus standard lead apron for radiation exposure of a simulated interventionalist, *JVIR* 2011; 22:437–442.
- Marx MV, Niklason L, Mauger EA. Occupational radiation exposure to interventional radiologists: a prospective study. *J Vasc Interv Radiol* 1992; 3:597–606.

36

## References (continued)

- Miller DL, Vano E, Bartal G, Balter S, Dixon R, Padovani R, Schueler B, Cardella JF, de Baere T. Occupational radiation protection in interventional radiology: A joint guideline of the Cardiovascular and Interventional Radiology Society of Europe and the Society of Interventional Radiology. *J Vasc Interv Radiol*. 21(5): p. 607-15.
- NCRP Report No. 116. Limitation of Exposure to Ionizing Radiation.. Bethesda, MD: National Council on Radiation Protection and Measurements, 1993.
- NCRP Report No. 122. Use of Personal Monitors to Estimate Effective Dose Equivalent and Effective Dose to Workers for External Exposure to Low-LET Radiation. Bethesda, MD: National Council on Radiation Protection and Measurements, 1995.
- NCRP Report No. 168. Radiation Dose Management for Fluoroscopically-guided Interventional Medical Procedures. Bethesda, MD: National Council on Radiation Protection and Measurements, 2011.
- Sanchez R, Vano E, Fernandez JM. Staff radiation doses in a real-time display inside the angiography room. *Cardiovasc Intervent Radiol* 2010; 33:1210-1214.
- Wörgul BV, Kundiyew YI, Sergiyenko NM, et al. Cataracts among Chernobyl clean-up workers: Implications regarding permissible eye exposures, *Radiat. Res.* 2007; 167(2), 233–243.
- World Health Organization. Efficacy and radiation safety in interventional radiology. Geneva: World Health Organization, 2000.